


# Reticulum

(w/o Uranyl nitrat)



1 = Placenta; 2 = Kidney; 3 = Skin; 4 = Pancreas

The Reticulum Kit is an often used staining in histopathology laboratories. It belongs to the so called special stains. The Kit reagents are standardized and a time saving procedure. It is finished in 20 minutes.

The Reticulum Kit is used for a specific staining and the detection of reticular fibres, little or premature fibres of collagens which contains carbohydrate groups. The staining based on the method of Snook and is used in the differential diagnosis certain tumor types.

ProTaq's Reticulum doesn't contain Uranyl nitrat.

ProTaq's Reticulum 401901200 100 Tests

