


6. Microbiology

Microbiology is the study of microscopic organisms and implies numerous sub-disciplines including virology, mycology, parasitology, and bacteriology.



6.1. Stains for Microbiology

Cat.-No.	Product Name	Format	Application	Quantity
401901799	ProTaq [®] Acridin Stain Kit (Acid-Fast-Bacteria)	Ready-to-use	Microbiology	2 x 125 ml
401901500	ProTaq [®] Iodine Solution	Ready-to-use	Microbiology	120 ml
401901900	ProTaq [®] Kinyoun Acid Fast	Ready-to-use	Microbiology	100 Tests
401901400	ProTaq [®] Modified Trichrome (Microsporidia)	Ready-to-use	Microbiology	2 x 250 ml

6.2. FITC Stain Kits for Microbiology

ProTaq[®] FITC Stain Kits for Microbiology contain

- ▶ Primary Antibody
- ▶ Secondary Antibody (Two Step)
- ▶ FITC Protein Blocking Agent
- ▶ FITC Mounting Medium

Cat.-No.	Product Name	Source	Format	Application	Quantity
300600200	ProTaq [®] Stain FITC Borrelia Burgdorferi	Rabbit	Ready-to-use	Immunofluorescence	100 Tests
300600900	ProTaq [®] Stain FITC Candida Albicans	Rabbit	Ready-to-use	Immunofluorescence	100 Tests
300630400	ProTaq [®] Stain FITC Chlamydia Pneumoniae	Mouse	Ready-to-use	Immunofluorescence	100 Tests
300600500	ProTaq [®] Stain FITC Chlamydia Trachomatis EBS	Goat	Ready-to-use	Immunofluorescence	100 Tests
300600600	ProTaq [®] Stain FITC Chlamydia Trachomatis EBS	Rabbit	Ready-to-use	Immunofluorescence	100 Tests
300600700	ProTaq [®] Stain FITC Chlamydia Trachomatis MOMP	Goat	Ready-to-use	Immunofluorescence	100 Tests
300630500	ProTaq [®] Stain FITC Giardia Lamblia	Mouse	Ready-to-use	Immunofluorescence	100 Tests
300600800	ProTaq [®] Stain FITC Hepatitis B Surface Antigen	Goat	Ready-to-use	Immunofluorescence	100 Tests
300600400	ProTaq [®] Stain FITC Legionella Pneumophila SG2	Mouse	Ready-to-use	Immunofluorescence	100 Tests
300600300	ProTaq [®] Stain FITC Leptospira Biflexa	Rabbit	Ready-to-use	Immunofluorescence	100 Tests
300630600	ProTaq [®] Stain FITC Pneumocystis Carinii	Mouse	Ready-to-use	Immunofluorescence	100 Tests
300601000	ProTaq [®] Stain FITC Varicella Zoster	Mouse	Ready-to-use	Immunofluorescence	100 Tests

6.3. Mounting Media

Cat.-No.	Product Name	Format	Application	Quantity
401901691	ProTaq [®] FungiStainMount Blue	Ready-to-use	Microbiology	15 ml
401901991	ProTaq [®] FungiStainMount Red	Ready-to-use	Microbiology	15 ml

6.4. Buffer

ProTaq[®] Digest Myco Bacteria and **ProTaq[®] Tb Buffers** are a good choice in mycobacterial diagnostic. Containing an alkali solution and N-acetyl-l-cysteine solution ProTaq[®] Digest Myco Bacteria is a special system for digestion of sputum in mycobacterial procedures. ProTaq[®] Tb Buffers are pre-sterilized buffers used for washing sputums after digestion and decontamination.

Cat.-No.	Product Name	Synonym	Format	pH	Application	Quantity
400400199	ProTaq [®] Digest Myco Bacteria		Ready-to-use		Microbiology	10 x 100 ml
402002099	Tb Buffer w. Tween 80	TRIS Borate Buffer (w/o. EDTA)	Ready-to-use	6.4	Microbiology	2 x 500 ml
402001899	Tb Buffer	TRIS Borate Buffer (w/o. EDTA)	Ready-to-use	6.8	Microbiology	2 x 500 ml